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**Direttiva 93/42/CEE, Allegato V, attuata in Italia con Dlgs. 46 del 1997/02/24**  
93/42/EEC Directive, Annex V, transposed in Italy by Dlgs. 46 of 1997/02/24

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This Certificate is subject to CERMET regulations and it is valid only for the above mentioned Medical Devices that are subject to survey. The technical sheet is an integrating part of this Certificate

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
**ATTENTION:** The operators must carefully read and completely understand the present manual before using the product.

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
#### Disposable resuscitator bag


The resuscitator bag Gima is conceived for use in addition to artificial respiration and cardiopulmonary reanimation.

#### Instructions for use

- Connect the oxygen tube to an adjustable oxygen source.
- In case it is connected to a reservoir, adjust the oxygen emission so that the reservoir completely expands during inhalation and almost completely empties during expiration thus recharging the bag.
- Before using the product, check for its efficiency preferably connecting it to a test lung: check that the connections, the reservoir, if present, and the valves complete the whole ventilation cycle correctly.
- In case of need connect the face mask to the apposite connector.
- Follow the necessary ventilation and cardiac resuscitation.
- Press the bag and observe the thorax expand during inhalation.
- Release the pressure on the bag and observe the thorax reduce during expiration.
- While ventilating check for:
  - Cyanosis symptoms
  - The correctness of the applied ventilation method
  - The airway pressure
  - The correct operation of all valves
  - The correct operation of the reserve and the oxygen tube

During ventilation the non-return valve could contaminate with vomit, blood or other secretions. In this case, stop ventilation and proceed as follows:

- Rapidly compress the bag several times so that to expel the substance through strong jets of air.
- If the non return valve is not clean yet: dip the valve in the water and rapidly compress the bag several times with strong jets of air. If the valve is still obstructed, throw out the bag.

### Prescriptions

- Do not use the bag in presence of toxic substances.
- Remove the oxygen reserve and the supplementary valve if oxygen is not being given. The non-removal reduces the recharge speed and consequently the ventilation frequency.
- Do not give supplementary oxygen in presence of free flames.
- Do not use oil, grease or other hydrocarbon containing substances on the bag and its components. The supplementary oxygen, given under pressure, can cause explosions if mixed with hydrocarbons.
- The use of this device is reserved to qualified medical personnel trained in lung ventilation and cardiac resuscitation advanced techniques.
- Before using the device on patients, make sure you are able to assemble, disassemble and use it.
- During ventilation, always check the airways pressure with a pressure gauge
- Qualified personnel trained in PEEP (Positive End Expiratory Pressure) only can use PEEP technique with this kind of device.

### Technical specifications

Bag capacity: 1600 ml adult mod. (34245-77); 500 ml child mod. (34246-8); 280 ml infant mod. (34247-9)  
Ejection capacity: 700 ml (one hand) – 900 ml (two hands) adult mod. - 300 ml - 350 ml child mod.  
- 150 ml - 225 ml infant mod.

Reserve capacity (optional): 2500 ml for adult and child models - 500 ml infant mod.

Ventilation cycle rate: 20 cycles/min. adult and child models - 40 cycles/min. infant mod.

Oxygen concentration: with reserve 99% - without reserve 45% for adult and child; 90% infant

Exerted pressure intervals:  $60 \pm 10$  cmH<sub>2</sub>O max. for adult mod.  
 $40 \pm 5$  cmH<sub>2</sub>O for child and infant mod.

The results have been obtained in the following conditions:

Adult: Vt-600ml, conformity 0.02L/cmH<sub>2</sub>O, resistance 20cm H<sub>2</sub>O/L/s

Child: Vt- 300, conformity 0.01L/cmH<sub>2</sub>O, resistance 20cm H<sub>2</sub>O/L/s

Infant: Vt-70ml, conformity 0.01L/cmH<sub>2</sub>O, resistance 20cm H<sub>2</sub>O/L/s

The bag performance can vary according to the operator and different factors such as room temperature, patient's lung capacity, ventilation rate and operator's hands size.

### Insufflateur de reanimation monousage

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