GIMA

PINZE BIPOLARI
BIPOLAR FORCEPS
PINCES BIPOLAIRES
PINZAS BIPOLARES
PINÇAS BIPOLARES
ΔΙΠΟΛΙΚΕΣ ΛΑΒΙΔΕΣ

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
MANUAL DE USO Y MANTENIMIENTO
MANUAL DE USO E MANUTENÇÃO
ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ

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FEATURES

The bipolar forceps GIMA are surgical instruments to be used with equipment for bipolar electrosurgery. All models are covered with insulating material, can be sterilized in autoclave and are equipped with universal coupling for use with Diatermo GIMA and the most common models on the market.

PRESCRIPTIONS

The product must be used by qualified medical personnel only, trained in electrosurgery techniques and maintenance operations on this kind of products.

Particular attention must be paid as regards the instruments care. The incorrect use or the lack of proper maintenance can make the instruments quickly deteriorate.

*These instruments can cause shocks or burns if used improperly!
Avoid putting these instruments in contact with others, which are not insulated!*

Do not use electro-surgery instruments on patients with Pacemakers, or in presence of flammable or anaesthetic liquids!

Diathermic knife electrodes used with these instruments cause the destruction of tissues and are extremely dangerous if they are not correctly used. It is therefore recommended to carefully read the operating instructions supplied by the diathermic knife electrode manufacturer.

The product must be sterilized before using it. It is recommended to sterilize it in autoclave.

UNPACKING

Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.

They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

MAINTENANCE

It is important to carefully examine the product before using it to monitor the presence of breaks, defects or wear and tear. Particular attention must be paid to moving parts, tips and connectors. The instruments must be kept in suitable boxes or protected trays, with their tips covered with gauzes or apposite caps in dry and clean place.

In case the product is defective or excessively worn out, do not use it and do not try to repair it.

For repair apply to your own retailer, GIMA or qualified personnel.

After using it, the product must be thoroughly cleaned with warm water, disinfected and rinsed possibly with distilled water, carefully dried and sterilized before the next use.

The product can be cleaned in apposite ultrasound tanks. Follow the instruction of ultrasound tank manufacturer.

The product can be disinfected in surgical instrument washing machine at a temperature up to 95°C (203°F).

After cleaning and before sterilizing them, it is advisable to oil the instruments with physiologically safe oil, in particular tips, connectors, terminals and all moving parts.


Be careful that the product does not get in contact with acids or aggressive disinfectants, which can corrode it.

The product is reusable and in conformity with ANSI/AAMI standards on sterilizing.

We guarantee the product for a minimum of 20 sterilizing cycles, if they are sterilized in compliance with the following criteria. Most instruments fully pass this number of cycles in presence of a care and maintenance protocol.

The recommended sterilizing method is steam sterilization in autoclave; the cycle temperature must not exceed 134°C (275°F) for up to 10 minutes, to avoid damaging the product.

It is recommended to follow the instructions of the autoclave manufacturer and arrange the instruments in such a way that their surface is directly and fully exposed to steam. Make them cool before using it. Hot air sterilization is recommended.

GUARANTEE

Thank you for having acquired our product. This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.

During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, etc. No compensation will be given for loss of use of the product.

Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.

The guarantee will be invalidated if the registration number has been removed, cancelled or altered. Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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