GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy

Italia: tel. 199 400 401 - fax 199 400 403

Export tel. +39 02 953854209/221/225 fax +39 02 95380056

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

Saturimetro Oxy-100 **Oxy-100 Pulse Oximeter**

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1. Safety

1.1 Instructions for Safe Operations

- Check the device to make sure that there is no visible damage that may affect user's safety and measurement performance. When there is obvious damage, stop using the device.
- Necessary service must be performed by qualified service engineers ONLY. Users are not permitted to repair it by themselves.
- The oximeter cannot be used together with the devices not specified in User Manual.

1.2 Cautions


- Explosive hazard—DO NOT use the oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- DO NOT use the oximeter while the testee is under MRI or CT scanning.
- The uncomfortable or painful feeling may appear if using the sensor of this device continuously on the same place for a long time, especially for the patients with poor microcirculation. It is recommended that the sensor should not be applied to the same location for over 2 hours. If any abnormal condition is found, please change the position of sensor.
- For the individual patients, there should be a more prudent inspecting in the placing process. The sensor can not be clipped on the edema and tender tissue.
- To dispose the device or its accessories, the local law should be followed.

1.3 Attentions

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- If the oximeter gets wet, please stop operating it. When it is carried from cold environment to warm and humid environment, please do not use it immediately.
- DO NOT operate the button on front panel with sharp materials.
- High temperature or high pressure steam disinfection to the oximeter is not permitted. Refer to related chapter for instructions of cleaning and disinfection.
- The intended use of this device is not for therapy purpose.

2 Overview

2.1 Appearance

Display screen: display SpO2 plethysmogram and parameter values.

Navigation key:

Up/Left/Increase

Press this key, the default screen can be shifted to display pulse rate (PR) or perfusion index (PI). When on the system setup screen, press it to move the cursor upwards or to the left and adjust parameter values.

Down/Right/Decrease

Its function is similar with the key: "Up/Left/Increase".

Mode/OK: press this key, the screen can be shifted between default screen and other screen layout; longtime press this key, the menu screen will be displayed; when you finish parameter setting, press this key to confirm.

Data interface: used for uploading data.

Reserved port: for future use.

Power/Back: Power on/off the device by longtime pressing; short time press it to back to upper level operation.

Recall/Backlight: Longtime press it to enter SpO2 recall screen; short time press it to turn on/off the backlight "SpO2": SpO2 Probe Connector.

2.2 Structure

It consists of the main unit and SpO2 probe.

2.3 Features

- It is light, small in size and easy to carry.
- High resolution LCD to display plethysmogram and measured data.
- Monitor SpO2 and PR parameters simultaneously;
- PI (Perfusion Index) display is available;
- Real-time store and recall SpO2 and Pulse Rate value.
- Audible & Visual alarm function;
- Power management function;

2.4 Scope of Application

Oxy-100 is applicable for monitoring SpO2 and PR.

2.5 Working Environment

Operating temperature: 5~40jæ

Operating humidity: 30~80%

Atmospheric pressure: 70kPa~106kPa

2.6 Installation of Battery and Holder

1. Open the rear panel with coin or an ordinary flat screwdriver, as shown in Figure 1.


Figure 1

2. According to the polarity mark, insert three AA batteries into battery house, as shown in Figure 2.


Figure 2

3. Close the battery cover and fasten it.

4. Fixing Holder (figure 3).


Figure 3

3. SpO₂ Probe Connection

Connect the SpO₂ probe to the right panel's connector labeled "SpO₂". After starting the monitor insert one finger into the probe (index finger, middle finger or ring finger with proper nail length) according to demonstration, as shown in the following figure.


Figure 5

Instructions of Operation

1. The finger should be put properly and correctly.
2. Do not shake the finger and keep at ease during using.
3. Do not put wet finger directly into sensor.
4. Avoid placing the sensor on the same limb which is wrapped with a cuff for blood pressure measurement or during venous infusion.
5. Do not let anything block the emitting light from the sensor.
6. Vigorous exercise and electrosurgical device interference may affect the measuring accuracy.
7. Using enamel or other makeup on the nail may affect the measuring accuracy.

4. Operation

4.1 Default Screen

Press “” power key for 2 seconds to start the Oximeter, then insert finger into rubber cushions of the probe, the screen will display the default screen, as shown in Figure.

Figure 6

Screen Description:

“%SpO2” “99”: SpO2 value; “99” is a percentage value;

“10:12:45”: Current time;

“PR 65”: Pulse rate value: 65bpm.

“”: Battery indication;

“”: Data memory icon; when the data on screen becomes stable, the device will start to store data automatically. No icon means the current data will not be stored. If the memory is full or the total number of the records is 256 pieces, the earliest records will be overwritten and the icon “” will appear on the screen for prompt, as shown in the figure.

Note: It is suggested that the data shall be uploaded to computer for saving, or the earliest records will be overwritten.

“”: Alarm Icon

“”: Pulse strength bar-graph

Lower part displays plethysmogram.

4.2 Display Screen with PI Value

On the default screen, press ▲ / ▼ Navigation key to shift screens between default screen and display screen with PI value. The display screen with PI value is shown below.


Figure 7

4.3 Other Screen Layouts

On default screen, press “” Mode/OK key to enter into other screen layouts (as shown in Figure 8A).

▲ / ▼ **(Navigation key)**: shift screens between PR observation screen (Figure 8A) and PI observation screen (Figure 8B).

 **(Mode/OK key)**: shift screens among four screens (as shown in Figure 6, Figure 7, Figure 8A and Figure 8B.)


Figure 8A PR Observation Screen


Figure 8B PI Observation Screen

4.4 Indication for No Signal

If “No Signal” prompts on screen (as shown in Figure 9), please check whether the probe connects well and whether the finger is in appropriate position.

Figure 9 Probe Off Screen

4.5 Menu Setup

On the above mentioned screens, longtime press “ ” key for entering into menu screen (as shown in Figure 10).


Figure 10 Menu Screen

Screen Description:

“**SETUP**”: set parameter values, refer to 4.5.1 for details.

“**UPLOAD DATA**”: enter into uploading status, refer to 4.5.2 for details.

4.5.1 SETUP

On the menu screen, select “**SETUP**” and then press “” key for entering into system setup screen. The setup screen is as shown in the following figures.

SETUP (Ver1.0)	
DATE	Oct. 18,07
TIME	10:12:45
CONTRAST	4

Figure 11 - System Setup Screen (A)

ALARM SOUND	ON
ALARM-SpO ₂	90%
ALARM-PR HI	120
ALARM-PR LO	50

Figure 11 - System Setup Screen (B)

ALARM-PR HI	120
ALARM-PR LO	50
PULSE SOUND	ON
REC INTERVAL	1s

Figure 11 - System Setup Screen (C)

Operation Instructions:

DATE: Date setting

1. When cursor stays on the Month of the date, press “” (Mode/OK) key to active Month option, the cursor flashes on the Month of the date;
 2. Press  /  (Navigation key) to adjust month.
 3. Press “” (Power/back) key or “” (Mode/OK) key to confirm and exit from date setting.
 4. The procedures of adjusting Day value and year value are the same with Month adjustment.
- Date Format: MM. dd, yy.

Note: *The setting operations of other parameters (such as TIME, CONTRAST, ALARM SOUND etc.) are the same with date setting.*

TIME: Time setting

CONTRAST: LCD contrast setting

ALARM SOUND: The factory default value is “ON”. When the set is “ON”, the device stays in alarm status and the alarm icon displays on the upper right display screen.

ALARM SpO₂: SpO₂ alarm setting; the factory default value: 90%.

ALARM PR HI/LO: Pulse rate: High/Low limit setting; the factory default value: High limit: 120, Low limit: 50.

PULSE SOUND: Pulse beep sound switch.

REC INTERVAL: Record interval setting; five options: “1s, 2s, 4s, 8s” and “OFF”; When the option is “OFF”, the device will not record real-time measurement data.

Note: *The length of data record is limited to 30 seconds at least, and the maximal length for one record is also limited to one hour (for one second interval) , 2hours(for 2 seconds interval) 4 hours (for 4 seconds interval) or 8 hours (for 8 seconds interval).*

4.5.2 UPLOAD DATA

On the menu screen, select "UPLOAD DATA" and then press " " key for entering into connecting status (as shown in Figure 12). When you transmit data (SpO2 and PR values) to your computer, please let the oximeter stay in connecting status.


Figure 12 Connecting Status Screen

4.6 Data Recall

On the default display screen, longtime press Recall/backlight key to enter into recall list screen.

Oct. 18,07	12:09:35
Oct. 18,07	15:07:35
Oct. 18,07	10:03:35
Oct. 18,07	12:50:35

Figure 13 Recall List

4.6.1 Data Recall

Choose one record in the recall list, then press " " key, the display screen will display trend graph, as shown in Figure 14A.


Figure 14A Trend Graph Recall Screen

Screen description:

“**SpO₂**”: The left ordinate value is %SpO₂ value;

“**12:50:35**”: Measuring time;

“**2s**”: Record interval is 2 seconds.

“**PR**”: The right ordinate value is pulse rate value.

Operation Instructions:

- press “” key to shift trend recall screens (as shown in Figure 14A £→Figure 14B and Figure 14C).


Figure 14B - Trend Graph Recall Screen


Figure 14C -Trend Graph Recall Screen

- Press “▲/▼” key to view trend graph.
- Press “⏻” key to return to recall list screen.

4.6.2 Data Deletion

Choose one record in the recall list, then longtime press “⚙️” key, an icon “🗑️” will be displayed in the back of this record, as shown in Figure 15.

Oct. 18,07	12:09:35	🗑️
Oct. 18,07	15:07:35	🗑️
Oct. 18,07	10:03:35	🗑️
Oct. 18,07	12:50:35	🗑️

Figure 15

When finishing the selection of the records needed to be deleted, press “⏻” key, then the screen will prompt that “Delete all the marked records?” as shown in Figure 16, at this time select “Yes” and press “OK” key to delete all the selected records.


Figure 16

5 Technical Specifications

A. Display mode:

dot-matrix monochrome LCD;

B. Power supply requirement:

Supply voltage: 3.6V~5.0V (3 xAA size) alkaline battery

Operating current: <50mA

Continuous working time: >30 hours

C. SpO2 Parameter Specifications

Measuring range: 35~99%

Accuracy: $\pm 2\%$ (during 75~99%)

$\pm 3\%$ (during 50~74%)

Default alarm limit: 90%

D. Pulse Rate Parameter Specifications

Measuring range: 30bpm~240bpm

Accuracy: ± 2 bpm or $\pm 2\%$ (whichever is greater)

Default alarm limits: High: 120bpm Low: 50bpm

E. Perfusion Index Display

Range: 0~20%

F. Data Record

Record data once every 1/2/4/8seconds, up to 70-hour records.

G. Resistance to interference of surrounding light:

The difference between the value measured in the condition of indoor natural light and that of darkroom is less than $\pm 1\%$.

H. Resistance to interference of man-made light:

Values of SpO2 and Pulse Rate can be accurately measured by pulse oxygen simulator.

I. Dimensions: 145 mm (L) x 74 mm (W) x 29 mm (H)

Net Weight: 210g (including batteries)

J. Classification

The type of protection against electric shock:

Internally powered equipment

The degree of protection against electric shock:

Type BF applied part.

The degree of protection against harmful ingress of liquids: Ordinary equipment without protection against ingress of water. Electro-Magnetic Compatibility:

Group I, Class B

6 Accessories

1. A probe
2. A holder
3. Battery AA x 3
4. A User Manual
5. A Quality Certificate
6. A data cable (optional)
7. Oximeter Data Manager software (optional)

Note: The accessories are subject to change. See the Packing List for detailed items and quantity.

7 Repair and Maintenance

7.1 Maintenance

The life of this device is 5 years. In order to ensure its long service life, please pay attention to the use of maintenance. A. Please change the batteries when the low-voltage indicator appears.

B. Please take out the batteries if the oximeter will not be used for a long time.

C. The best storage environment of the device:

Ambient temperature: -20 to 60°C

Relative humidity: 10~95%

7.2 Cleaning and Disinfecting Instruction

A. Surface-clean sensor with a soft gauze by saturating with a solution such as 75% isopropyl alcohol.

B. Then Surface-clean with clean water and dry with a piece of cloth.

*High-pressure sterilization cannot be used on the device.
Do not immerse the device in liquid.*

***Disposal:** The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of

correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

8 Troubleshooting

Trouble	Possible Reason	Solution
The SpO2 and Pulse Rate display instable	<ol style="list-style-type: none"> 1. The finger is not placed inside enough. 2. The finger is shaking or the patient is moving. 	<ol style="list-style-type: none"> 1. Place the finger properly and try again. 2. Let the patient keep calm
Can not turn on the device	<ol style="list-style-type: none"> 1. The batteries are drained or almost drained. 2. The batteries are not inserted properly. 3. The device's malfunction. 	<ol style="list-style-type: none"> 1. Change batteries. 2. Reinstall batteries. 3. Please contact the local service center.
No Display	<ol style="list-style-type: none"> 1. The device will power off automatically when there is no signal and no operation for 1 minute. 2. The batteries are almost drained. 	<ol style="list-style-type: none"> 1. Normal. 2. Change batteries.

Appendix

Symbols on the screen

Symbol Description

%SpO₂ The pulse oxygen saturation

PI Perfusion Index

♥ bpm Pulse rate (Unit: beats per minute)

 Low battery voltage

 Alarm Icon

 Memory Icon

Symbols on the panels

Symbol Description

 With Type BF applied part

SpO₂ SpO₂ probe connector

 Power/Back Key

 Mode/OK Key

 Recall/Backlight Key

 Navigation Key

 Data Interface

SN Serial Number

 Separate collection for this device