**ARTICOLI E APPARECCHI PER MEDICINA
PROFESSIONAL MEDICAL PRODUCTS**

**ASPIRATORE TOBI
ASPIRATEUR TOBI
TOBI ASPIRATOR
ASPIRADOR TOBI
ABSAUGER TOBI**

**MANUALE D'USO
MODE D'EMPLOI
INSTRUCTION MANUAL
MANUAL DE INSTRUCCIONES
HANDBUCH**

GIMA Spa Via Trento, 10 - 20064 Gorgonzola (MI) Italia
Tel.02 9517251-2-3-4 Fax.02 95300275 E-Mail gimaita@tin.it
Export Tel. ++39 02 9517251-2-3-4 Fax++39 02 95304706 E-Mail gimaexp@tin.it

ISTRUZIONI

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
INSTRUCTIONS

- TOBI surgical aspirator with 230V~ power supply
- Easily transportable from one hospital ward to another, particularly suitable for tracheotomized patients, small surgical applications and post-operative therapy at home.
- To be used for aspirating body fluids such as, for example, mucus or catarrh.
- Easily portable equipment designed for virtually continuous use.

Made of highly heat-resistant, electrically insulated plastic material in conformity with the latest European safety standards.

Supplied with a polycarbonate 1-litre autoclavable vessel complete with overflow valve.

Equipped with aspiration regulator, located on the front panel.

PRECAUTIONS

WARNING

**CAREFULLY CONSULT THE INSTRUCTION MANUAL
BEFORE USING THE EQUIPMENT**

WARNING

THE EQUIPMENT MAY ONLY BE UTILIZED BY QUALIFIED PERSONNEL

WARNING

**DO NOT OPEN THE EQUIPMENT. FOR TECHNICAL ASSISTANCE
CONTACT GIMA AFTER-SALES SERVICE**

GIMA does everything within its power to ensure that its products are safe to use.

It is, however, of fundamental importance to observe the basic rules of safety in order to avoid electric shocks or fires.

1. On using the equipment for the first time, and on each successive occasion, check that the power cord cable is not damaged, split or peeling. **If it is, do not insert the plug in the electric socket**
2. On using the equipment for the first time, check that no parts of it have been broken, cracked or opened during transportation. **If so, do not insert the plug in the socket**
3. Do not handle the plug with wet or damp hands. Do not use the equipment in the shower or in the bath.
4. Do not immerse the equipment in water
5. Do not allow children or incapable people to use the equipment without supervision
6. Do not leave the equipment exposed to atmospheric agents or heat sources
7. **Do not utilize the equipment for uses other than those foreseen by the manufacturer**

Technical characteristics

Model	TOBI
Dimensions	370mm X 220mm X 210mm
Weight	2.5 Kg
Power supply	230V 50 Hz
Absorbed power	230VA
Maximum suction	-0.70 bar
Maximum suction flow rate	15 L/min
Fuse	2 x 1.6 A

SYMBOLS

	Class 2 isolation equipment
	CE marking in conformity with EC directive 93/42/CEE
	Warning, consult the instruction manual
	Type B equipment
	Fuse
230V ~	Alternate current
50 Hz	Mains frequency
IP33	Degree of protection against ingress of liquids

ACCESSORIES SUPPLIED

DESCRIPTION
COMPLETE VESSEL 1000 ml
SINGLE VESSEL 1000 ml
VESSEL COVER
SEALS FOR 1000 AND 2000 ml VESSELS
CONICAL FITTING
TUBING SET
EQUIPMENT TO VESSEL CONNECTING TUBE
TUBE WITH PROBE
ASPIRATION PROBE
ANTI-BACTERIAL FILTER

CLEANING OF ACCESSORIES

Wearing disposable latex gloves, proceed to clean the accessories as follows:

1. Grip the vessel with your hand and turn the cap in a counter-clockwise direction.
2. Remove the red cage covering the float located in the vessel cap (conical part with 3 openings)
3. Once the cage is removed, the float with the red seal comes out of the vessel cap.
4. **To disinfect the accessories, use either denatured alcohol or hypochlorate-based solutions, easily purchasable from chemist shops.**
5. Once cleaning has been completed, insert the float in the cage, with the seal facing towards the cage opening.
6. Insert the opening of the cage into the seat in the vessel cap.
7. To facilitate the cleaning of the vessel, fill it 1/3 full of ordinary water. The aspirated liquid, thus diluted, will then be easier to remove.

THE ACCESSORIES MAY BE DISINFECTED USING DENATURED ALCOHOL OR HYPOCHLORATE-BASED SOLUTIONS

PERIODICAL MAINTENANCE CHECKS

THE TOBI EQUIPMENT HAS BEEN DESIGNED TO LIMIT MAINTENANCE OPERATIONS TO THE MINIMUM. IT SHOULD HOWEVER BE CHECKED DAILY PRIOR TO USE.

1. Hook up the power cord to the electricity mains
2. Switch the machine on
3. Close the aspirator outlet with your finger
4. Check, to the touch, that the equipment is aspirating properly.
5. Rotate the knob from right to left and check, to the touch, that the aspiration regulating control is in perfect working order.


N°	Type of fault	Cause	Remedy
1	No aspiration	Vessel cap badly screwed down	Unscrew the cap, then rescrew it correctly
2	No aspiration	Cap seal not in its seat	Unscrew the cap and insert the seal properly in its seat
3	No aspiration	Folds in external tubes	Arrange the tubes so as to prevent kinking.
4	The float does not close	If the cap has been washed, ensure that the float is not partially detached	Snap the float into place
5	Slow aspiration		Fill the vessel to 1/3 full of ordinary water
6	No aspiration due to outflowing of mucus	Filter blocked	Replace filter
	Faults 1 - 2 - 3 - 4 - 5	None of the remedies has achieved the desired results	Contact the seller or GIMA After-sales Assistance Service

IN THE EVENT OF A FAILURE OR MALFUNCTIONING DO NOT ATTEMPT TO REPAIR THE EQUIPMENT, BUT CONTACT GIMA AFTER-SALES SERVICE IMMEDIATELY.

IF, ON INSPECTION OF THE EQUIPMENT BY GIMA AFTER-SALES SERVICE, THE EQUIPMENT IS FOUND TO HAVE BEEN TAMPERED WITH, THE GUARANTEE WILL BE CONSIDERED NULL AND VOID.

OPERATING INSTRUCTIONS

- Connect the short silicone tube with the filter 1 to the aspirator outlet 2; the other end of the tube must be connected to the outlet with the vessel cap float 3.
- Connect the long silicone tube 3 to the other outlet of the vessel cover (the outlet without the float).
- Attach the probe connection ⑤ to the other end of the tube 5
- Connect the plug to the electrical mains supply
- Switch on the equipment 6

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