MOPEDIA

MA RP5XX 01 A_02-2014

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I Class Medical Device according to directive 93/42/EEC
and subsequent amendments

1. Codes

- RP520 4-wheel aluminum rollator with seat ATLANTE
- RP525 4-wheel aluminum rollator with seat and footrest GAYA
- RP530 4-wheel aluminum rollator with seat OCEANO
- RP535 4-wheel aluminum rollator with seat URANO

2. Introduction

Thank you for purchasing the rollator, **MOPEDIA** by MORETTI S.p.A., realized to help people with movement difficulty inward and outward, during voyages or in everyday life. The production has been made according to the highest quality and safety standards in accordance with European Directive on Medical devices (93/42/EEC).

This instruction manual gives you suggestion for a correct use and for your safety. We recommend you to read carefully this instruction manual before using. In case of doubts, please contact your dealer.

Note: Check all parts to verify it won't be any shipping damage. In case of shipping damages, do not use. Contact the dealer for further instructions.

3. Intended use

This rollator fit for people with movement problems. Can be used as assistance to facilitate ambulation or even servo-assisted ambulatory. Anyway, it can't be used by people with evident malfunction of the hand or the arm, even by people with balance problems. The rollator can be used inward and outward. We recommend using on flat surface where the wheels can reach the land. It can't be used on rock ground or on slanted area as the rollator can be unstable. The device is not intended to carry people or things.

WARNING!

Do not use the product for a purpose not indicated in this manual. Moretti Spa declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product. The manufacturer can modify the device or the instruction manual without prior notification.

4. Declaration of conformity EC

Moretti S.p.A. declare under sole responsibility that the product made and traded by Moretti S.p.A. and belonging to rollator range are in accordance with the European Directive on Medical devices (93/42/EEC). Moretti S.p.A. guarantees and declared under sole responsibility what follow:

1. The devices under discussion satisfy the essential requisites requested by the I° enclosure directive 93/42/CEE as laid down by the VII° enclosure of the above mentioned Directive.
2. The devices under discussion ARE NOT MEASURING INSTRUMENTS.
3. The devices under discussion ARE NOT MADE FOR CLINICAL PROBES.
4. The devices under discussion are traded with NON STERILE BOX.
5. The devices under discussion belong to class I°.
6. IT IS NOT ALLOWED to use and install devices except for the use provided by Moretti S.p.A.

Moretti S.p.A. places technical documentation at relevant authorities in order to prove the conformity of 92/42/EEC device.

7. Moretti S.p.A maintains and provides for authorities, the technical documentation proving the accordance to European Directive on Medical devices (93/42/EEC).

5. Warnings

5.1 General warnings

- Do NOT install or use the device before reading and understanding this manual; If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional dealer or technical personnel before attempting to install this equipment to avoid injure or damage;
- The rollator should only be used by patients with a suitable motor capacity which ensures their use without risks;
- Advise with physician or therapist to determinate the correct adjustment and the correct use;
- Physician or therapist should contribute to height adjustment for the maximum support, for a correct brake activation and for a general assembly;
- The rollators are not intended to be used as self-propelled devices.

5.2 Warnings for use

- Insure the brake cable to be correctly set;
- Verify the correct operation of the brake after every height adjustment to avoid a free flexion of the cable;
- Breaking this simple rules can cause personal damages;
- The rollator can be used only on a hard, flat surface. Do not use the rollator to climb the stairs
- Leaning always on both dumbbells to maintain balance;
- Pay particular attention to wet, icy or snowy;
- Do not overload the rollator;
- Do not use the rollator seat (where present) to carry persons or property, the seat is for the user's session with the parking brake engaged;
- All wheels must be in contact with the floor at any time during use to ensure the proper balancing of the rollator and maximum safety;
- The brake must be engaged before using the seat.

Note When the rollator is used in steady position, the handbrake must be locked

5.3 Stability warnings

- The rollator can provide ambulatory assistance to an individual weighing up to 136 kg, INCLUDING the loading bag;
- The rollator bag has a weight limitation of 5 kg.

5.4 Assembling & folding warnings

- Verify that both handle support tubes are correctly inserted in the frame and secured with locking knobs locked into the holder with its main wing nuts and knobs without movement;
- Both handle tubes are securely locked in the main support part with wing nuts and knobs without any movements;
- Do not hanging anything on the rollator frame to avoid distortions of the tubes when the rollator is open and damage the stability. Injury may also be caused by distortions;
- Take care to entrapment of fingers and / or clothing during the opening and closing.

6. Before use

Before using the device, always check:

1. That all parts of the rollator are locked in a secure manner.
2. Check that all screws are tight.
3. Make sure that the seat is properly locked
(if that 'not press further with the power seat bottom frame).
4. Check that the locking of the seat is correctly inserted.
5. Check the correct operation of the brake and the regular brake.

7. Assembling instruction

7.1 Before assembling the product

Unpack the product and check that no damage has occurred during shipping. In the event when the rollator is damaged, please contact your nearest dealer.

7.2 Assembling model RP520

- Open the rollator pulling the front part by leveraging on the legs.
- Place the handle tubes in the frame.
- Turn the locking knobs until they are locked.
- The hexagonal hole is needed to accommodate the screw head

7.2.1 Folding the rollator

Make lever on the lateral bars and bring back the front of the rollator.

7.2.2 Adjustment of the handles

Adjust the height of the handles using the side knobs. Unscrew the knob, place the handles to the desired height, turn the knob clockwise to tighten and fasten it again. Repeat for the other handle and make sure that both handles are positioned at the same height.

WARNING!

Make sure that the handle is in line with the direction of movement of the rollator. If the handle does not stay straight, there is a risk that the rollator wobbling or tipping over, especially if you push it harder

7.2.3 Basket Storage

The shopping basket is hung by two hooks to the structure of the rollator. (only model RP520)
The basket is placed in the lower part of the structure of the rollator.

WARNING! Use the basket for small items only. Do not overload the basket to prevent imbalances and the risk of overturning

7.2.4 Seat

Before sitting, place the rollator on a flat and solid surface and then lock the brakes. Make sure the latch of the lock is secured properly and that the rollator is in a stable position.

WARNING! Questions to sit down with your weight on the seat assembly. Do not sit by pushing only on one side

7.3 Assembling model RP525

- The rollator is dispatched with footrest packed separately.
- Grasp the two knobs and exert a force toward the outside so as to distend the session. (the first time that the rollator will be mounted handles with the knobs will not be connected, and then to open the rollator grasp the two tubes the uprights).
- Press the seat of the rollator up to arrest the frame tubes in position lock and then lock the safety clip
- Before use, always make sure that the rollator does not tend to close up under the effect of a small lateral force.
- Install the footrest on the rollator frame

WARNING! Make sure that the knobs is located into the holes

7.3.1 Folding the rollator

Lift up the seat firmly pulling the straps in the middle of seat.

7.3.2 Installation of the handles and adjustable handles

Insert the handles into the frame. Simply depress the spring button and set to the desired height. Repeat the opposite side and make sure that both handles are secure and set at the same height before use. The rollator is equipped with height adjustable handles through the pins to meet all user needs.

WARNING!

Adjust the height of the rollator so that a user who is standing with his arms at his sides have the pulse height handles. This will ensure the positioning of the arms bent at about 20° -30° when using the rollator

7.3.3 Bag installation

To install the bag, hang the ears of the bag on the seat, stop the tape with Velcro to the structure

WARNING! Use the basket for small items only. Do not overload the bag to prevent imbalances and the risk of overturning

7.3.4 Seat

Before sitting , place the rollator on a flat and solid surface and then lock the brakes. Make sure the latch of the lock is secured properly and that the rollator is in a stable position.

WARNING! Questions to sit down with your weight on the seat assembly. Do not sit by pushing only on one side.

7.4 Assembling model RP530

- Turn the wheel in the correct direction, unbuckle the belt which locks the rear wheels and lock the rear wheels and front with knobs supplied
- Grasp the two knobs and exert a force toward the outside so as to distend the session. (the first time that the rollator will be mounted handles with the knobs will not be

- connected, and then to open the rollator grasp the two tubes the uprights).
- Press the seat of the rollator up to arrest the frame tubes in position lock and then lock the safety clip
- Before use, always make sure that the rollator does not tend to close up under the effect of a small lateral force.

WARNING! Make sure that the knobs is located into the holes

- Before use, always make sure that the rollator does not tend to close up under the influence a small lateral force

WARNING! Be careful not to twist the brake cables when connecting the handlebar and wheels

7.4.1 Folding the rollator

Lift up the seat firmly pulling the straps in the middle of seat

7.4.2 Installation of the handles and adjustable handles

Insert the handles into the frame. Simply depress the spring button and set to the desired height. Repeat the opposite side and make sure that both handles are secure and set at the same height before use. The rollator is equipped with height adjustable handles through the pins to meet all user needs.

WARNING!

Adjust the height of the rollator so that a user who is standing with his arms at his sides have the pulse height handles. This will ensure the positioning of the arms bent at about 20 ° -30 ° when using the rollator.

7.4.3 Bag installation

To install the bag, hang the ears of the bag on the seat, stop the tape with Velcro to the structure.

WARNING! Maximum capacity for bag 5 kg. Use the bag for small items only. Do not overload the bag to prevent imbalances and the risk of overturning

7.4.4 Seat

Before sitting , place the rollator on a flat and solid surface and then lock the brakes. Make sure the latch of the lock is secured properly and that the rollator is in a stable position.

WARNING! Questions to sit down with your weight on the seat assembly. Do not sit by pushing only on one side.

7.5 Assembling model RP535

- Open rollator
- Release the red safety lock, if it is inserted.
- Grasp the two knobs and exert a force toward the outside so as to distend the session. (the first time that the rollator will be mounted handles with the knobs will not be connected,

- and then to open the rollator grasp the two tubes of the vertical uprights).
- Press the seat of the rollator up to arrest the frame tubes in position lock and then lock the safety clip.
- Insert the rear wheels in the structure of the rollator by means of snap button.

WARNING!

Beware the button is outside the structure to ensure proper placement of the wheel

- Before use, always make sure that the rollator does not tend to close up under the effect of a small lateral force.

WARNING!

Be careful not to twist the brake cables when connecting the handlebars and wheels

- Place the front wheels in the frame, the wheels are properly seated when you hear a “click”. Remove the wheels by pressing the red button on the structure of the rollator.

7.5.1 Fold the rollator

- Lift up the seat firmly to close the rollator
- Insert the red safety lock when closed completely the rollator

7.5.2 Installation of the handles and adjustable handles

Insert the handles into the frame, and lock them in an appropriate manner. The rollator is equipped with height adjustable handles through the pins to meet all user needs.

WARNING!

Adjust the height of the rollator so that a user who is standing with his arms at his sides have the pulse height handles. This will ensure the positioning of the arms bent at about 20 ° -30 ° when using the rollator.

7.5.3 Installation of the backrest

Insert the backrest into the guides on the plastic structure of the rollator

7.5.4 Bag installation

To install the bag, hang the ears of the bag on the seat , stop the tape with Velcro to the structure

WARNING! Maximum capacity for bag 5 kg. Use the bag for small items only. Do not overload the bag to prevent imbalances and the risk of overturning

7.5.5 Seat

Before sitting , place the rollator on a flat and solid surface and then lock the brakes. Make sure the latch of the lock is secured properly and that the rollator is in a stable position.

WARNING! Questions to sit down with your weight on the seat assembly. Do not sit by pushing only on one side.

8. How to use

- Start pushing when you are in an upright position with both hands on the handle;
- Make sure your hands are dry before using the rollator to avoid the risk of slipping;
- Walking maintaining an upright position with his eyes facing forward and not down;
- Do not push the rollator too far away from the body, this may cause falls and loads inadequate;
- Do not use the rollator to overcome steps higher than 20mm;
- When using outdoors, prefer the sidewalks properly prepared in accordance regulations for the removal of architectural barriers;
- To overcome small obstacles to leverage the rear footrests to raise the front of the rollator (models RP525 - RP530 - RP535)

9. Brakes

9.1 To brake

To brake the rollator pull the brake lever with your fingers upwards until the rollator slows down or stops. To continue to move the rollator, release the brake handles.

9.2 Lock the brake (handbrake)

Press down on the brake lever with the palm of your hand, when you hear a “click” wheels are locked. Lift the handle of the brake releasing it from the stop position.

WARNING! When the user is stationary, the brakes must be locked

10. Care and maintenance

The devices **MOPEDIA** by Moretti SpA on the market are carefully checked and carry the CE mark. For the safety of the patient and the doctor is recommended to be checked by the manufacturer or by an approved laboratory, at least every 2 years, the suitability for use of your rollator. If repairs are to be used only original spare parts and accessories.

The useful life of the product depends on the use, maintenance, storage and cleaning.

11. Cleaning

To clean parts, use only a damp cloth with a neutral detergent, not use cleaners containing solvents and metal brushes or abrasive or sharp objects, may damage the product.

12. Wheel care

1. Clean the wheels with warm water and mild detergent. Wipe dry with a clean cloth.
2. Periodically lubricate the wheels.

13. Technical features

Series ATLANTE										
Product codes	RP520									
Seat	38x34 cm									
Frame	Aluminum									
Weight and dimensions										
RP520	69	84 99	56	61	44	36	7,4	20	20	136
Series GAYA										
Product codes	RP525									
Seat	45x31 cm									
Frame	Aluminum									
Weight and dimensions										
RP525	50	86 98	56	67	44,5	27	9,7	20	20	136
Series OCEANO										
Product codes	RP530									
Seat	45x20 cm									
Frame	Aluminum									
Weight and dimensions										
RP530	69	79 92	53,5	62	46,5	27x74h	7,6	20	17,5	136
Series URANO										
Product codes	RP535									
Seat	44x24 cm									
Frame	Aluminum									
Weight and dimensions										
RP535	72	86 96	55	65	48	29	7,2	20	17,5	136

14. Spare parts / accessories

For spare parts and accessories please refer to the Moretti S.p.A. master catalogue

15. Warranty

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be apply in the event damages caused by improper use, abuse or alteration, and for the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual.

Moretti is not responsible of consequent damages, personal injuries or whatever caused by or in relation with mistaken installation or with device improper use.

Moretti does not guarantee for damages caused by the following condition: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, alteration not authorized, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual.

The warranty doesn't cover components subject to wear and tear during the correct use of the device.

15.1 Repairing

Warranty repair: If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will estimate with customer if the fault can be cover from warranty. Moretti, at its unquestionable discretion, could repair or replace the item, with a Moretti dealer or with own base. Labor cost can be in charge to Moretti if the repair is cover from warranty. A repair or a replacement doesn't extend the warranty.

Repair not cover from warranty: A product without warranty can be given back after Moretti authorization. The labor and shipping costs for a bath lift out of warranty are to be payed by the customer or by the dealer. The repairs are guaranteed for 6 months as from the good received.

Non-defective device: The customer will be informed if, after the device return and examination, Moretti declares that the device is not faulty. In this case the good will be sent back to customer, the shipping charge will be to be payed by the customer.

15.2 Spare parts

The original Moretti spare parts are guaranteed for 6 months as from the receiving date.

15.3 Exempt clauses

Moretti doesn't offer any other declarations, warranty or conditions, explicit or implicit, including possible declarations, warranties or condition of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti doesn't guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be fixed by according to the law is limited by the warranty period, in the limits of law. Some states or countries don't allow to limit the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user.

The present warranty may be modified without prior notification.

WARRANTY CERTIFICATE

Product _____

Purchasing date _____

Authorized dealer _____

Street _____ **Place** _____

Sold to _____

Street _____ **Place** _____

 **MORETTI S.p.A.**

Via Bruxelles, 3 - Melegnano 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com

MADE IN P.R.C.

Via Bruxelles, 3 - Melegnano
52022 Cavriglia (Arezzo)
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200
www.morettispa.com
info@morettispa.com