MAG1000

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User manual

Mode d'emploi

I.A.C.E.R. Srl
Via S. Pertini 24/A – 30030
Martellago (VE) ITALY

Tel. +39 041 5401356 • Fax +39 041 5402684
e-mail: iacer@iacer.it <http://www.itechmedicaldivision.com>

Elettromagnetoterapia modello

MAG
1000

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Electromagnetotherapy model

MAG 1000

<i>Summary</i>	22
<i>Introduction</i>	23
<i>Electromagnetotherapy</i>	23
<i>Technical Specifications</i>	24
<i>Manufacturer</i>	24
<i>Declaration of Conformity</i>	24
<i>Specifications:</i>	25
<i>Purpose</i>	25
<i>Technical specifications</i>	26
<i>Labelling</i>	27
<i>Labels details</i>	27
<i>Symbols</i>	28
<i>Kit contents</i>	29
<i>How to use</i>	29
<i>Warnings</i>	29
<i>Electromagnetic interference</i>	30
<i>Side effects</i>	30
<i>Instructions</i>	31
<i>Stored programs</i>	33
<i>Therapeutic belt positioning</i>	34
<i>Battery charging</i>	34
<i>Replacing the battery</i>	35
<i>Maintenance</i>	35
<i>Cleaning</i>	35
<i>Carriage and storage</i>	35
<i>Disposal</i>	36
<i>Troubleshooting</i>	36
<i>Repairing</i>	37
<i>Spare parts</i>	37
<i>Warranty</i>	37

Electromagnetotherapy

Electromagnetotherapy distinguishes itself from the “traditional” magnetotherapy since it uses high frequency electromagnets instead of permanent magnets: particularly, we refer to electromagnetotherapy when using devices emitting high frequency and low intensity pulsed electromagnetic fields (PEMF, carrying frequency from 20 to 30 MHz; with frequencies up to 5.000 Hz).

The principal use is focused on pain therapy: its action is to restore a cellular biomagnetic field that the illness has weakened.

The use of high frequency and to low intensity PEMFs, allows to get ample therapeutic results without collateral effects. This also allows the use in acute pathological trials. The electromagnetotherapy is particularly suitable for the care of the soft tissues pathologies, with extraordinary results in the regeneration of the tissues themselves.

Due to its characteristics, the electromagnetotherapy is universally recognized as the most suitable technique for the treatment of the bony pathologies, in particular for the osteoporosis.

*There are lots of effects that can be brought back to electromagnetotherapy: the piezoelectric effect, the effect of orientation of the collagene, the stimulation of the calcic deposition (**Barker - Lunt 1983, Bassett – Pawluk – Pilla 1974, Bassett - Valdes – Hernandez 1982**).*

Since the present day, all the equipments of electromagnetotherapy were built particularly to satisfy or the medical sector (with high prices) or the economic demands of patients, but with products of low quality.

MAG1000 is born really to reconcile the demand to have an equivalent device, for performances and effectiveness, to those devoted to medical sector, also maintaining a simplicity of use and an extremely favorable price.

Manufacturer

I.A.C.E.R. S.r.l.

Via S. Pertini, 24/a • 30030 Martellago (VE) ITALY

Tel. 041.5401356 • Fax 041.5402684

LACER S.r.l. is an Italian medical devices manufacturer (CE medical certificate n° MED24021).

Declaration of Conformity

LACER S.r.l., headquartered in Italy, via S. Pertini 24/A 30030 Martellago (VE), declares on its own responsibility that MAG1000 is manufactured in conformity with Directive 93/42/EEC (MDD), Annex II as modified by Directive 2007/47/CE dated 5 September 2007 (D. Lgs. 37/2010 dated 25 January 2010).

Notified Body: Cermet, Via di Cadriano 23 – 40057 Cadriano di Granarolo (BO) Italy

MAG1000 is a Class IIa equipment, with reference to Directive 93/42/EEC (MDD), annex IX rule 9 (and following modifications).

Martellago, 21/03/2010

*Legal representative
Mario Caprara*

Specifications:

MAG1000 has the following specifications:

- *Class IIa equipment (Directive 93/42/EEC, annexed IX rule 9 and following modifications);*
- *Class II, applied part type BF (Classif. EN 60601-1);*
- *Equipment not protected against liquids penetration;*
- *Equipment and accessories not subjected to sterilization;*
- *Use of the equipment is prohibited close to flammable substances or in environments with high concentrations of oxygen;*
- *Continuous operating mode equipment;*
- *Equipment not suited to be used in external.*

Purpose

Clinical purpose: *Therapeutic*

MAG1000 pathologies applications:

- *wrist articulation*
- *hand articulation*
- *shoulder articulation*
- *foot articulation*
- *ankle articulation*
- *knee articulation*
- *skeletal motor apparatus*
- *arthrosis*
- *atrophies and muscular distrofia*
- *bursitis*
- *bruises*
- *degeneration of locomotor apparatus*

- *distorsions*
- *neuralgias*
- *periarthritis*
- *benign lesions and muscular tears*
- *tendonitis*

MAG1000 is particularly suitable for the treatment and the care of the osteoporosis and all the pathologies on bony tissue, for dermatologic applications (varicose ulcers, phlebitis ulcers, arteropatic ulcers, acne, burn sore, not healed wounds, cellulitis, cutaneous illnesses), as well as to accelerate the recovery of scalpel cuts, to regenerate tissues, to avoid that wounds leave unsightly scars.


MAG1000 is recommended for domestic use

Expected lifetime: 5 years

Technical specifications

<i>Power supply</i>	<i>Rechargeable batteries, Ni-MH AA1600mAh 4,8V</i>
<i>Battery charger</i>	<i>XKD-C0300NHS6.0-0.5W, input 110-240VAC, 50-60HZ; output 6VDC 300mA</i>
<i>Max current abs.</i>	<i>≤ 150 mA (therapy mode)</i>
<i>Insulation class (CEI EN 60601-1)</i>	<i>II</i>
<i>Applied part (CEI EN 60601-1)</i>	<i>BF</i>
<i>Dimensions (mm)</i>	<i>134×64×18</i>
<i>Field intensity</i>	<i>Adjustable on 3 levels (low, medium, high)</i>
<i>Squared wave frequency</i>	<i>from 8Hz to 640Hz</i>
<i>Environmental conditions of operation</i>	
<i>Temperature</i>	<i>from +5 to +40 °C</i>
<i>Humidity</i>	<i>from 30 to 80%</i>
<i>Pressure</i>	<i>from 500 to 1060 hPa</i>

Labelling

Labels details

Label 1

ELECTROMAGNETOTHERAPY

MODEL: MAG 1000 
 S/N: MA5011/1 000001

MAX INPUT: DC 4.8V 200mA

OUTPUT POWER: 2x40mW (1 Kohm Load)

INDUCTIVE POWER: 2x10mW (150 ohm Load)

FREQc: 27 MHz


FREQm: 8-640Hz

I.A.C.E.R. Srl, Via S. Pertini 24/A
 30030 Martellago (VE) – ITALY

     
 0476 2008-11


Compliance with MDD 93/42/EEC

Label 2


Charger XKD-C0300NHS6.0-5W


Label 3

Label 4

Symbols

	Operating instructions
	Product subject to WEEE regulations concerning separate waste collection of electronic equipment
	Class II equipment
	Applied part type BF
	Non Ionizing radiation
	Compliance with Directive 93/42/EEC (MDD) and following modifications
	Manufacturing date

Kit contents

n° 1 MAG1000 device

n° 1 battery charger;

n° 1 elastic therapeutic belt;

n° 1 MAG1000 user manual.

How to use

Warnings

MAG1000 does not require any technical or professional knowledge to be use, and so it can be used both by professionals (physiotherapists, medics) and home users.

Take care of what follows:

- *Take care of position and meaning of the labels on MAG1000;*
- *Do not damage the therapeutic belt acting on connection cable, and avoid to roll up the cable around the belt or the device;*
- *Check the integrity of the battery charger before use; avoid the use in case of damage to the case or to the wire;*
- *Avoid the use of MAG1000 to people not educated through the reading of the manual;*
- *Avoid the use of MAG1000 contemporary with ointments containing free ions of magnetic metals;*
- *Do not use MAG1000 in damp environments;*
- *Do not use in presence of inflammable agents.*
- *Do not wear metallic objects during therapy;*
- *Take care to place the green side of the therapeutic belt on the skin.*
- *Use only cables and applicators supplied by the manufacturer. Inadequate cables and applicators could damage the device and/or could be hazardous for the patient.*
- *The user must periodically verify cables and applicators insulation and control their integrity.*

WARNING. Connect the battery charger plug to 230V main only when connected to MAG1000 device for batteries recharging. Disconnect the battery charger from main after each use.

WARNING. During therapy it is possible to hear a light hiss coming from MAG1000: such operation is normal and must not arouse worry.

The manufacturer is considered responsible of the performances, reliability and safety of the instrument only if:

- *possible additions, changes and/or reparations are effected from authorized personnel;*
- *the electric plant of the environment in which MAG1000 is inserted it is conforming to the national laws;*
- *the instrument is employed in conformity to the instructions contained in this manual.*

Electromagnetic interference

MAG1000 doesn't produce and doesn't receive interferences from other equipments. However, it's opportune to use the instrument keeping it at least at 3 meters from televisions, monitors, cellular telephones or any other electronic equipment.

The functioning of other devices connected to the patient can be compromised by the functioning of short wave therapy devices.

Contraindications and side effects

Patient in pregnancy, tuberculosis, juvenile diabetes, viral (in acute phase) illnesses, mycosis, cardiopathic subjects, serious arrhythmias or pacemaker carriers, children, metallic prosthesis carriers, acute infections, epileptics (different medical prescriptions excepted).


The functioning of some electric implantable devices, such as pacemakers, could be compromised during a treatment with a short wave therapy device. Consult a physician before starting the treatment.

No significant side effects are known of, nor are reported particular contraindications for excessive time length using the device.

Instructions

- 1) Switch on MAG1000 pressing ON/OFF key  ;
- 2) Check battery status: if the icon  on the upper side of the display is flashing, proceed to recharge the unit following the instructions (see chapter "Battery charging");
- 3) Connect the therapeutic belt (or belts) to the plug A/B on the side of MAG1000;
- 4) Chose therapy program running through programs menu, using PR+ and PR- keys;
- 5) Set therapy time using TM+ to increase and TM- to decrease: each program has its own set therapy time, which can be modified increasing or decreasing it with 5 minutes steps;
- 6) MAG1000 is now ready for use: the display shows informations about the program in use, the working frequency (not adjustable), the therapy time, the treatment area (or areas, indicated by icons near the human body on the display) and field intensity;
- 7) Press  key : a sound signal will advise the therapy starting, the display will show the flashing icons near human body and the remaining therapy time, while it will be possible to modify field intensity during the whole treatment, pressing more times IN. key. The relative icon is positioned on the lower right side of the display, with the indication of the three intensities (low intensity-one bar, medium intensity-two bars and high intensity-three bars);

- 8) The therapy can be suspended pressing the key   : the flashing icons of treatment areas will become fixed and the display will show a flashing icon  on the upper right side. To restart therapy press again   ;
- 9) To end the therapy press key  : MAG1000 will shown again the selected program and setted therapy time. To switch off MAG1000, keeping pressed key  ;
- 10) A sound signal inform of the end of the program, and the display will show the flashing icons "0'00"; to start again a new therapy , press key  and then follow again the indications from step 4);

WARNING: the display shows two icons relative to channel A and B. The flashing icon shows the not in use channel (no therapy belt connected).


Stored programs

Pr. N°	Freq. (MHz)	ON/OFF Modulation (Hz)	Time (h/min)	Name	Area
P.1	27	8	45 min.	Arthritis-Arthrosis	Body, arms, legs
P.2	27	8	2 h	Rheumatisms	Body, arms, legs
P.3	27	8	4 h	Pain- Articular prosthesis	Body, arms, legs
P.4	27	16	45 min.	Herniated disc	Body
P.5	27	16	2 h	Sciatica	Body, legs
P.6	27	16	4 h	Osteoporosis	Legs
P.7	27	32	2 h	Fractures	Body, arms, legs
P.8	27	32	4 h	Sprains	Arms, legs
P.9	27	48	2 h	Dislocations	Arms, legs
P.10	27	48	4 h	Traumas / Bruises	Arms, legs
P.11	27	64	2 h	Myalgia (tears/contractures/cramps)	Arms, legs
P.12	27	64	4 h	Tendonitis	Arms, legs
P.13	27	80	2 h	Hematoma	Body, arms, legs
P.14	27	80	4 h	Microcirculatory insufficiency	Body, arms, legs
P.15	27	160	2 h	Phlebitis	Arms, legs
P.16	27	160	4 h	Varicose veins/Varix Peripheral circulation	Arms, legs
P.17	27	320	2 h	Scars / Dermatitis / Dermatosi	Body, arms, legs
P.18	27	320	4 h	Allergies	Body, arms, legs
P.19	27	640	2 h	2h Medium application	Body, arms, legs
P.20	27	640	4 h	4h Long application	Body, arms, legs

IACER S.r.l. suggest the therapy time as shown on the table and pre-setted on MAG1000. However, the user can adjust the time as he prefers. MAG1000 uses therapy time values, working frequency values and field intensity values coming from scientific and medical literature, as result of well known sperimentations and clinical evaluations (Barker - Lunt 1983, Bassett - Pawluk - Pilla 1974, Bassett - Valdes - Hernandez 1982).

Therapeutic belt positioning

Here below a list of main positions for the therapeutic belt.

Battery charging

If the battery status is insufficient to complete the working program, the display will show the icon  on the upper left side, near wording PROGRAM.

To proceed with the charging, switch off MAG1000 keeping pressed  and connect the battery charger to the plug on the upper side of MAG1000.

A complete recharge of the batteries needs 8/10 hours.

MAG 1000 is equipped with a software protection that switch off the equipment when the battery charger is connected.

Replacing the battery

If, after normal recharging, the battery cannot sustain a full treatment session, recharge again, as it is probably approaching the end of its useful life. If even after a second recharge the battery still cannot complete a treatment session, follow the steps here below:

- *Contact the manufacture or authorized dealer for original spare part;*
- *Open the battery compartment on the back side of MAG1000, disconnect the red/black wire from connector and remove the exhausted battery;*
- *Insert new battery and connect the red/black wire to the plug;*
- *Close battery compartment*

WARNING: *for disposal exhausted battery follow the instructions on chapter “Disposal”. Do not open or burn battery. Do not short circuit poles. Keep battery away from sparks or naked flames. In the event of internal electrolyte coming into contact with skin or garments, wash immediately with water. In the event of electrolyte coming into contact with eyes, rinse thoroughly and seek medical attention.*

Maintenance

Cleaning

Clean the equipment from the dust using a dry soft cloth.

Resistant stains can be removed using a sponge soaked in solution of water and alcohol (20%).

Carriage and storage

Carriage precautions

MAG1000 is a portable device, so it does not need any particular carriage precautions.

However we recommend to put away MAG1000 and its accessories in their own bag after every treatment.

Storage precautions

MAG1000 is protected till following environmental conditions:

Outside of the packaging

<i>Temperature</i>	<i>from +5 to +40 °C</i>
<i>Humidity</i>	<i>form 30 to 80%</i>
<i>Pressure</i>	<i>form 500 to 1060 hPa</i>

Inside of the packaging

<i>Temperature</i>	<i>from -10 to +55 °C</i>
<i>Humidity</i>	<i>from 10 to 90%</i>
<i>Pressure</i>	<i>from 500 to 1060 hPa</i>

Disposal

The equipment is subjected to WEEE regulations (see the symbol  on the label) concerning separate waste collection: when disposing this product, please use the designed areas for disposing electronic waste or contact the manufacturer.

Troubleshooting

If it is used in accordance with the instructions of the user manual, MAG1000 does not need a particular regular maintenance.

If you find any malfunctioning using MAG1000, please follow these instructions:

- *Check battery status (see chapter “Battery charging”);*
- *Check the correct connection between MAG1000 and the therapeutic belt;*
- *Check all the operations have been done properly;*
- *Every two years we suggest a complete check of device (contact the manufacturer or locator dealer).*

WARNING: after prolonged inactivity, the unit may not switch on. This is because the battery is completely drained and not due to malfunction. Recharge the battery as described above before troubleshooting for faults.

Repairing

Every intervention on device must be performed by manufacturer or authorized personnel.

Spare parts

Contact the manufacturer for original spare parts.

To preserve product warranty, functionality and product safety we recommend to use only original spare parts.

Warranty (we recommend that you fill in and post the light blue warranty card).

MAG1000 is covered by a 24-month warranty with effect from the date of purchase. The batteries, battery charger, connection cables and therapeutic belt are covered by a 6-month warranty; parts subject to wear are excluded from the warranty, unless evident manufacturing defects are found.

The warranty will lapse in the event that the equipment is tampered with or work has been carried out on the same by personnel not authorised by the manufacturer or authorised dealer.

The warranty conditions are described in the “Warranty Regulations”.

Post the light blue card and keep the green card

N.B. responsibility of the purchaser: the purchaser is responsible for validating the warranty by filling in the warranty card (light blue warranty validation card), stamped by the retailer, and sending it to the manufacturer.

In the event of future repairs under warranty, the equipment must be packaged to prevent damage during transport and sent to the manufacturer together with all of the accessories, ensuring that all parts of the service card (green service card) have been filled in. The purchaser only has the right to repair under warranty when the equipment is returned to the manufacturer complete with the green service card signed by the retailer, together with the receipt or invoice proving the correct origin of the product.

Warranty regulations

1. *We recommend that the purchaser fills in the light blue warranty validation card, stamped by the retailer, and posts it to the manufacturer immediately after purchase.*
2. *In the event of repairs under warranty, the purchaser must fill in all parts of the green service card and send it to the manufacturer together with the equipment. The receipt for tax purposes or the purchase invoice must be sent together with the goods.*
3. *The electronic parts are covered by a 24-month warranty. The warranty is given through the point of sale or directly from the manufacturer.*
4. *The warranty covers exclusively product damage causing operational defects.*
5. *The warranty covers exclusively the repair or replacement free of charge, including labour, of components found to be defective in terms of manufacture or material.*
6. *The warranty does not apply to damage caused by neglect or use not complying with the instructions provided, damage caused by work carried out by unauthorised personnel or damage caused by accidental causes or the buyer's negligence, with particular reference to external parts.*
7. *The warranty does also not apply to damage to the equipment caused by incompatible power supplies (the equipment operates at 4.8V with an internal battery pack).*
8. *Parts subject to wear after use are excluded from the warranty.*
9. *The warranty does not include transport costs to be paid by the purchaser in relation to the method and speed of transport.*
10. *The warranty expires after 24 months. After such time repair work will be carried out at the rates currently in force for the parts replaced and the labour and transport costs.*
11. *Any controversy will fall within the exclusive jurisdiction of the Venice courts.*

MAG1000. All rights reserved. MAG1000 and the  logos are owned by I.A.C.E.R. Srl and are registered.

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MAG1000. All rights reserved. MAG1000 and the  logos are owned by I.A.C.E.R Srl and are registered.
I.A.C.E.R S.r.l. SU