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Venduto a _____

Via _____ Località _____

Hangzhou Sejoy Electronics & Instruments Co. Ltd Building 2, No.202, Zhenzhong Road,
West Lake Economy & Technology Zone 310030 Hangzhou China

EC REP

Shanghai International Holding Corp. Gmbh (Hamburg) - Eiffestrasse 80, 20537 Hamburg, Germany

Distrib. MORETTI S.P.A. Via Bruxelles, 3 - Meleto 52022 Cavriglia (Ar) Italy - www.morettispa.com

INDEX

1. CODE	pag.13
2. INTRODUCTION	pag.13
3. INTENDED USE.....	pag.13
4. DECLARATION OF CONFORMITY.....	pag.13
5. GENERAL WARNINGS.....	pag.13
6. ELECTROMAGNETIC COMPATIBILITY WARNINGS	pag.14
7. SIMBOLOGY	pag.14
8. GENERAL DESCRIPTION.....	pag.14
8.1 Parts view	pag.14
9. USE WARNINGS.....	pag.15
10. ADVANTAGES	pag.15
10.1 Body temperature.....	pag.15
11. HOW TO USE	pag.16
11.1 Temperature measurement	pag.16
11.2 Temperature taking hints	pag.16
11.3 Memory mode.....	pag.17
11.4 Replace the battery.....	pag.17
11.5 After measurement	pag.17
11.6 Selecting a temperature scale	pag. 17
12. MAINTENANCE.....	pag. 17
13. CLEANING AND DISINFECTION.....	pag.17
13.1 Cleaning	pag.17
13.2 Disinfection	pag.17
14. CONDITIONS OF DISPOSAL.....	pag.17
15. SPARE PARTS AND ACCESSORIES.....	pag.18
16. TROUBLESHOOTING.....	pag.19
17. ELECTROMAGNETIC COMPATIBILITY DECLARATION	pag.19
18. TECHNICAL FEATURES	pag.21
19. WARRANTY	pag.21
20. REPAIRING.....	pag.21
21. SPARE PARTS.....	pag.21
22. EXEMPT CLAUSES.....	pag.22

CE 0197 -IIa Class Medical Device

According to Directive 93/42/EEC and further modifications

1. CODE

RM350 Infrared ear thermometer

2. INTRODUCTION

Thank you for purchasing a GRADO by Moretti Thermometer. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

3. INTENDED USE

The thermometer is a device designed to measure body temperature via eardrum.

CAUTION!

- Do not use the product for a purpose not indicated in this manual;
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product;
- The manufacturer reserves the right to change the information contained in this document without previous notice.

4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under his sole responsibility that the product placed in the market by the same Moretti S.p.A. is in conformity to the relevant provisions of Directive 93/42/CEE on medical devices and subsequent amendments and are produced by Hangzhou Sejoy Electronics & Instruments Co. Ltd according to the standards indicated in the regulations:

- CEI 60601-1 - Part 1: General requirements for basic safety and essential performance
- EN 60601-1-2 - Part 1: General requirements for safety - Collateral standard:
Electromagnetic compatibility
- ISO 80601-2-56:2009:Medical electrical device - (Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for measuring body temperature).
- Manufacture are ISO 9001 and ISO 13485 certified.

5. GENERAL WARNINGS

- Read instructions thoroughly prior to use;
- Store and use thermometer away from mobile phones;
- Always take temperature readings from the same ear, as ears may differ in temperature;
- Keep the battery away from children;
- Keep the unit dry and away from any liquids and direct sunlight;
- Infrared thermometers should not be used shortly after exercise, bathing or coming indoors;
- Choking Hazard: Thermometer cap and battery may be fatal if swallowed. Do not allow children to use this device without parental supervision;
- Do not store or keep in direct sunlight, or in high temperature because the packaging is made in paper;
- Shelf life - the thermometer use limit is defined by the wear of the parts;
- Do not allow children to play on or operate the thermometer;
- Do not place thermometer battery near extreme heat as it may explode;
- It is recommended the performance should be checked every two years;
- Remove battery from the device when not in operation for a long time;

- The use of temperature readings for self-diagnosis may lead to the worsening of existing disease conditions;
- Do not attempt measurements when the thermometer is wet as inaccurate readings may result;
- Do not bite the thermometer. Doing so may lead to breakage and/or injury;
- Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings;
- After each use, disinfect the thermometer especially in case the device is used by more than one person;
- Disinfettare il termometro dopo ogni uso, specialmente nel caso di riutilizzo su più pazienti.

6. ELECTROMAGNETIC COMPATIBILITY WARNINGS

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2:2007. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:


- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Consult the manufacturer or field service technician for help.

CAUTION!

- Do not use in proximity of the device, mobile phones and other devices, which can generate strong electrical or electromagnetic fields as they can interfere with the operation of the device


7. SIMBOLOGY

8. GENERAL DESCRIPTION

The thermometer is specifically designed for safe use in the eardrum. The Infrared Ear Thermometer is a device capable of measuring the people's body temperature by detecting the intensity of infrared light emitted from the ear canal of human. It converts the measured heat into a temperature reading displayed on the LCD. The unusual backlight design is to the benefit of reading in the night.

8.1 Parts view

1. Probe
2. ON-OFF
3. Test
4. Battery Cover

Fig.1

9. USE WARNINGS

CAUTION!


-  The performance of the device may be degraded should one or more of the following occur:
- Operation outside the manufacturer's stated temperature and humidity range;
 - Storage outside the manufacturer's stated temperature and humidity range;
 - Patient temperature is below ambient temperature.

- Do not bend, drop the thermometer;
- Do not boil probe to disinfect;
- Do not disassemble, except in case of battery replacement;
- The thermometer should only be used under the supervision of an adult. Prohibit children from biting, hitting and removing the battery cover, dislodging and potentially swallowing the battery;
- Do not use the device in case of strong electromagnetic field such as microwave oven or cell phone.

10. ADVANTAGES

Thermometry measures core body temperature, which is the temperature of a body's vital organs. (Figure 2) Eardrum temperatures accurately reflect core body temperatures since the eardrum is in the cranial cavity and is not affected by the environment. The eardrum shares blood supply with the hypothalamus (part of the brain.) Therefore, body temperature changes are reflected sooner in the ear than they are in other parts of the body.

Vital Organs

- **Axillary:** (armpit) temperature measures skin temperature and therefore may not be a reliable indicator of core body temperature.
- **Oral:** temperature can be influenced by drinking, eating, and breathing.
- **Rectal:** temperature readings often lag behind changes in core body temperature. A cross-contamination risk may exist as well.

Fig.2

10.1 Body temperature

- Body temperature can vary from 35°C to 38°C. The generally accepted "normal" temperature is 37.0°C;
- Usually your temperature is lower in the morning than in the afternoon. Normal temperature varies according to the location of the temperature reading.

The following table demonstrates how normal temperature ranges can vary by location of the body. Therefore measurements from different body locations should not be directly compared.

Normal Ranges by Body Location:

Axillary	34.7 - 37.3°C	94.5 - 99.1°F
Oral	35.5 - 37.5°C	95.9 - 99.5°F
Ear	35.8 - 38.0°C	96.4 - 100.4°F

A person's normal temperature range tends to decrease with age.

Normal Ranges by Age:

0-2 Years	36.4 - 38.0°C	97.5 - 100.4°F
3-10 Years	36.1 - 37.8°C	97.0 - 100.0°F
11-65 Years	35.9 - 37.8°C	96.6 - 99.7°F
Over 65 Years	34.8 - 37.5°C	96.4 - 99.5°F

Normal Ranges vary from person to person and can fluctuate throughout the course of a day.

11. HOW TO USE

CAUTION!

Temperature readings are available in the Fahrenheit or Celsius scale ($^{\circ}\text{C}/^{\circ}\text{F}$); located in the upper right corner of LCD.) With the unit off, press and hold the "On/Off" button for approximately 3 seconds to change the current setting. Release the "On/Off" until Diagram 3 is present, the unit is now ready for testing in the chosen scale.

11.1 Temperature measurement

1. Press the Mem Button. The LCD screen will appear as shown (Figure 3). Wait for the screen to change to setting mode as shown (Figure 4) and the flashing Arrow mark appears, "Bi---Bi---" sound is heard. The thermometer is now ready for a new measurement. If there is a mistake, the "Er1" or "Er2" appears;

Fig. 3


Fig.4

2. Place the probe snugly into the ear canal. Make sure the probe tip faces the eardrum as depicted (Figure 5). Press the Test Button. A long beep will be sounded 1 second after the "Test" Button is activated. This signals the end of the process as results are displayed on the LCD screen. Meanwhile, the backlight function will operate to help with reading in the night;
 3. The thermometer is ready for testing again once an arrow mark displays on the LCD and "Bi---Bi---" sound is heard. A waiting period of 20 seconds between testing is recommended to ensure complete accuracy;
4. The thermometer will shut off automatically after 60 seconds of inactivity. To prolong battery life, press and hold the Mem Button approximately 3 seconds to turn the unit off;


Fig.5

11.2 Temperature taking hints

1. Temperatures of two eardrums may slightly differ. Always test the same ear drum when comparing results.
2. External factors may influence ear temperatures such as
 - Lying on one ear or the other
 - Exposed to very hot or cold temperatures
 - Recent swimming or bathing

NOTE: Under these circumstances, remain in the testing environment for 20 minutes or longer prior to temperature taking.

3. The ear must be free from obstructions such as excess earwax;

Fig.6

Fig.7

4. When using the thermometer on infants under age 1, pull the ear up making sure the sensor faces the eardrum; (Figure 6)
5. When using the thermometer on individuals over the age of 1, pull the ear back making sure the sensor faces the eardrum; (Figure 7)
6. Do not touch the probe tip. A cotton swab and rubbing alcohol may be used to clean the sensor tip;
7. If thermometer is stored in a significantly different environment

- than testing location, place it in the testing location for approximately 30 minutes prior to use;
- Keep the ear canal clean.

N.B: If the temperature is below 32 ° C, the display will show "Lo", if the temperature exceeds 42.9 ° C, the display will show "Hi". In case the temperature will be over 37.8 ° C, the acoustic signal will have a very fast intermittently.

11.3 Memory mode

- The Memory Mode can be accessed two ways
 - After the thermometer has been turned on and displays the screen depicted in Diagram 3, press the Mem Button, the mark "O" will appear in the upper left corner of the display. (Figure 8)
 - After testing, press the Mem Button the mark "O" will appear in the upper left corner of the LCD.
- The thermometer will automatically memorize the last 10 temperature readings. Each time the Mem Button is pressed, the screen displays past readings that correspond with a number 1-0. The number 1 reflects the most recent reading while the number 0 reveals the oldest reading stored in memory.
- In the memory mode, arrow mark always exists. The user can press the Test Button to take new measurements.

11.4 Replace the battery

- Battery needs replacing when the LCD displays the icon shown (Figure 8);
- Slide battery cover down as shown (Figure 9);

Fig.8

Fig.9

Fig.10

- Remove battery as shown (Figure 10);
- Replace with one new battery (CR 2032). Place the "+" side upward;
- Put back the battery cover;
- Take the battery out in case of long-term non-use.

11.5 After measurement

- Power off: Device will automatically shut off if left idle for more than 1 minute to extend battery life.
- Clean the probe after each use to ensure an accurate reading and avoid cross contamination. (See the section of Care and Cleaning for details.)

11.6 Selecting a temperature scale

- Temperature readings are available in the Celsius (°C) or Fahrenheit (°F) scale.
- With the unit off, press and hold the TEST BUTTON for 3 seconds to enter into unit changing mode.
- Press and release TEST BUTTON to select the unit.
- When the preferred unit on the display, press  BUTTON to exit the unit changing mode.

12. MAINTENANCE

The GRADO devices by Moretti S.p.A are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair, use only original spare parts and accessories.

13. CLEANING AND DISINFECTION

13.1 Cleaning

- The probe tip is the most delicate part of the thermometer. To safely clean the tip, gently wipe its surface with a cotton swab or soft cloth moistened with rubbing alcohol. Allow the probe to fully dry for at least 3 minutes;
- Use a soft, dry cloth to clean the thermometer display and exterior;

- The thermometer is not waterproof do not submerge the unit in water when cleaning;
- Put in the box and store thermometer in a dry location free from dust and direct sunlight.

13.2 Disinfection

If you need to disinfect the device use isopropyl rubbing alcohol.

NOTE: Never use acids, alkalis or solvents such as acetone

14. CONDITIONS OF DISPOSAL

General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the thermometer in the appropriate disposal areas for recycling.

 **WARNINGS FOR THE PROPER DISPOSAL OF THE PRODUCT, AS REQUIRED BY THE EUROPEAN DIRECTIVE 2002/96/EC:** At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin the product bears the barred bin symbol.

 **DISPOSAL OF WASTE BATTERIES (DIRECTIVE 2006/66/EC):** This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

15. SPARE PARTS AND ACCESSORIES

For spare parts and accessories refer to main catalogue.

16. TROUBLESHOOTING

Situation	Solution	Error message
Button malfunctions	Remove old battery and replace with new one	No
Temperature results are not within test range 32.0 °C~43.0°C (89.6 °F~109.4 °F) Lo = too low; H = too high	Make sure the probe tip and lens filter are clean. Be sure thermometer is properly inserted into a clean ear canal. Contact seller if error message continues after several consecutive attempts.	Lo H
Ambient temperature is not within the allowed operating range 10° C	Place thermometer in a room for 30 minutes with a temperature range between	Er 1
Other error	Contact seller	Er 2

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EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF Energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR11	Class B	The device is suitable for use in all establishments.
Harmonic emissions IEC61000-3-2	N/A	Including domestic establishments and those directly connected to the public low-voltage power supply network
Voltage fluctuations/flicker emissions IEC61000-3-3	N/A	

Guidance and manufacturer's declaration - Electromagnetic immunity -

The devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	± 2kV for power supply line ± 1kV for input/output line	N/A	
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	N/A	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	N/A	
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.

Note: UT is the a.c. mains voltage prior to the application of the test level

Guidance and Manufacturer's Declaration - Electromagnetic Immunity -

The devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM banda	N/A	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5G MHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey c, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 3 Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device

NOTE 4 Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 KHz to a 80MHz $d=1.2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

The user can prevent electromagnetic interference by maintaining a minimum distance between the thermometer and the device that might causes interference.

18. TECHNICAL FEATURES

Measuring site:	Ear canal (ear mode)
Measuring range:	32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
Accuracy:	± 0.2°C (0.4°F) during 35.5°C~42.0°C(95.9°F~107.6°F) at 16°C~35°C (60.8°F~95.0°F) operating temperature range ±0.3°C (0.5°F) for other measuring and operating temperature range
Display:	Liquid crystal display (lcd) with backlight 0.1°C or 0.1°F
Response time:	2 seconds
Operating temperature range:	10.0°C~40.0°C (50°F ~ 104°F); ≤85% RH, non-condensing
Storage and transportation temperature range:	-25°C to 55°C (-13°F to 131°F) ; ≤95% RH, non-condensing
Clinical repeatability:	±0.3°C(0.5°F)
Dimension (LxWxH):	118x 38x 35 mm
Weight (including battery):	Approx gr40(with batteries)
Battery:	One CR2032 lithium cell battery
Battery life:	Approx 1 year/3000 readings

19. WARRANTY

Every Moretti product is guaranteed for 2 (two) YEARS from the purchase date, except the restrictions specified below. The warranty shall not apply in case of damage caused by improper use, abuse, alterations of the product or any other use not in compliance with the user manual. The proper use of the product is indicated in the user manual.

Moretti S.p.A. shall not be liable for any damages or personal injuries due to incorrect assembly/improper use of the device, not compliant with the user manual. Moretti does not guarantee its products against damages or faults in the following cases: natural disasters, non authorized maintenance or reparations, damages caused by electrical supply problems (if provided), use of parts or components not provided by Moretti, failure in accomplishing the guidelines and instructions for use, non authorized alterations, damages due to the transport (different from Moretti original delivery) or failure in maintaining the device as indicated in the manual.

The components subject to wear are not covered by this warranty if the damage is due to the normal use of the product.

20.REPAIRING

Warranty repair

If a Moretti product shows material or manufacturing faults during the warranty period, Moretti S.p.A., together with the client, will consider whether the product fault is covered by warranty. Moretti S.p.A., at its unquestionable discretion, may replace or repair the item under warranty in a Moretti reseller point or at its own qualified offices. Labour costs for repairing the product may be at Moretti's expense if it is determined the reparation is under warranty. A reparation or replacement does not renew nor extend the warranty.

Repair not covered by warranty

A product, not covered by warranty, may be returned for reparation only after receiving prior authorization from Moretti client service. Labour and delivery costs for an out-of-warranty reparation will be fully at the client's or reseller's expense. Reparations of out-of-warranty products are guaranteed for 6 (six) months, from the day the client receives the repaired product.

Non defective device

The client is informed if, after examining and testing a returning product, Moretti decides the product is not defective. The product will be sent back to the client and all the return delivery costs will be at his expense.

21. SPARE PARTS

The original Moretti spare parts are guaranteed for 6 (six) months from the day you receive them.

22. EXEMPT CLAUSES

Except what clearly specified in this warranty and in accordance with the law, Moretti S.p.A. does not offer its clients any other declaration, warranty or condition, expressed or implicit, including any declaration, warranty or condition of merchantability, of non violation or non interference, of suitability for a specific aim. Moretti S.p.A. does not guarantee that the use of Moretti products will be continuous and faultless. The duration of possible implicit warranties that may be applied in accordance with the law is limited to the warranty period. Some states or countries do not allow limitations to the duration of an implicit warranty, or the exoneration, or the limitations of accidental or indirect damages related to products for consumers. In such states or countries some exoneration or limitation of this warranty may not be applied to the user. This warranty is subject to modifications without any prior notice.

**WARRANTY CERTIFICATE**

Product _____

Purchased on (date) _____

Retailer _____

Address _____ **Town/city** _____

Sold to _____

Address _____ **Town/city** _____

Hangzhou Sejoy Electronics & Instruments Co. Ltd Building 2, No.202, Zhenzhong Road,
West Lake Economy & Technology Zone 310030 Hangzhou China

EC

REP

Shanghai International Holding Corp. GmbH (Hamburg) - Eiffestrasse 80, 20537 Hamburg, Germany

Distrib. MORETTI S.P.A. Via Bruxelles, 3 - Meleto 52022 Cavriglia (Ar) Italy - www.morettispa.com

MORETTI S.P.A.

Via Bruxelles, 3 - Meleto
52022 Cavriglia (Arezzo)

Tel. +39 055 96 21 11
Fax. +39 055 96 21 200

www.morettispa.com
info@morettispa.com