MANUALE D'USO

SEDIA DOCCIA BASCULANTE

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USER MANUAL TILTING SHOWER CHAIR ACQUA OCEAN TRANSIT

INTRODUCTION

Thank you for choosing Demarta Virginio. The products have been designed and manufactured to meet all your needs for practical, correct and safe use. This manual contains small suggestions for the correct use of the device you have chosen and valuable advice for your safety. We recommend that you read this manual carefully before using the item you purchased. If in doubt, please contact the retailer, who will be able to help you and advise you correctly. The manufacturer reserves the right to make changes to the device and the following manual without notice in order to improve its characteristics.

GENERAL WARNINGS

Use the product in the manner indicated in this manual. If you still have doubts, contact DEMARTA VIRGINIO S.r.l.

Always evaluate the compatibility of the product with the characteristics of the subject together with your doctor or therapist.

Always check the integrity of the product before each use, even when new. Do not use if damaged or worn.

DEMARTA VIRGINIO S.r.l. declines all responsibility for damage to property or persons resulting from improper use, not agreed with a professional, not corresponding to what is indicated in the user manual.

Any accidents, anomalies, breakages must be reported to DEMARTA VIRGINIO S.r.l. in order to monitor the quality and critical issues of the product.

Keep the packaged product away from any heat source as the packaging is made of cardboard.

Always pay attention to the presence of children.

The product has been designed only for the intended use. Do not abuse it in any way.

EC DECLARATION OF CONFORMITY

DEMARTA VIRGINIO S.r.l. declares under its sole responsibility that the product manufactured and placed on the market by DEMARTA VIRGINIO S.r.l. complies with the applicable provisions of EU Directive 2017/745 on medical devices.

To this end, DEMARTA VIRGINIO S.r.l. guarantees and declares under its sole responsibility the following:

The devices meet the general safety and performance requirements as required by Annex I of Regulation 2017/745 as prescribed by Annex IV of the aforementioned regulation.

- The devices ARE NOT MEASURING INSTRUMENTS.
- The devices ARE NOT INTENDED FOR CLINICAL INVESTIGATIONS.
- The devices are marketed in NON-STERILE packaging.
- The devices are to be considered as belonging to class I in compliance with the provisions of Annex VIII of the aforementioned regulation.

INTENDED USE

Evacuation and mobility aid intended for people with walking difficulties in the bathroom.

SPECIFIC WARNINGS

Always be very aware of moving parts that could cause limb entrapment and personal injury.

Always consult your doctor or therapist for proper use of the device.

The service life of the device is determined by wear and tear of non-repairable and/or replaceable parts.

TECHNICAL FEATURES

Length (cm)	Width (cm)	Height (cm)	Seat Width (cm)	Seat Depth (cm)	Seat Height from Floor (cm)	Wheel Dimensions (cm)	Maximum Load (kg)	Net Weight (kg)	Box Dimensions (cm)
									
85	48	106	46	41	61 / 65,5	12.5 (5")	120	17	67 x 61 x 66

ASSEMBLY

1) Wheel assembly.

Remove the packaging frame and insert the 4 wheels into the appropriate holes, fixing them with the 4 plastic clips supplied.

⚠ ATTENTION!

After fitting the 4 wheels, make sure that they are properly stable before proceeding with the Mounting the frame.

2) Seat assembly.

Remove the seat from the packaging.

Place it on top of the frame and apply light pressure on it until the 4 clips are correctly locked.

ATTENTION!

Make sure it is properly secured before proceeding with assembly.

3) Inserting the seat end cap.

Remove the seat cap from the packaging and insert it into the hole provided between the basin and the seat.

4) Tray assembly.

Remove the basin from the packaging and place it under the seat.

Slide it in from the back of the chair, taking care to insert it correctly into the rails.

5) Backrest assembly.

Remove the backrest from the box.

Insert the two side tubes of the backrest inside the two ends of the frame.

Insert the star locking knobs into the holes provided on the sides of the frame, correctly fixing the backrest to the frame.

⚠ ATTENTION!

Screw the 2 star locking knobs all the way in, ensuring proper frame stability.

6) Headrest mounting.

Remove the headrest from the packaging.

Insert the headrest tube into the holder already installed on the push handle and secure it with the screw already on the frame.

ATTENTION!

The headrest can be easily adjusted to different positions. To make adjustments, unscrew one at a time the knobs located on the headrest and adjust their position. Then tighten the knob is good. Perform the same operations on the other knob if necessary.

7) Chest pads assembly.

Remove the chest pads from the packaging.

Attach the chest pads to the 2 ends of the frame, using the 4 Phillips screws provided.

ATTENTION!

The pair of lateral thrusts have different holes, which allow you to adapt to the different sizes of the shower chairs.

8) Footpeg assembly.

Remove the two footrests of the shower chair from the packaging.

The pull-out footpegs should be inserted into the two ends of the chair frame by inserting the black pin into the hole in the frame tube.

After inserting them correctly, push down on the footplate and rotate it towards the inside of the chair until it locks on the stop point.

To remove the footrest, use the release lever then rotate the footrest outwards and pull it upwards.

BEFORE EACH USE

- Always check the state of wear and tear of the mechanical parts involved so as to ensure that the device can be used in total safety for people and property.
- Check that the device is assembled correctly.

WARNINGS FOR USE

- Do not exceed the maximum permitted capacity (120 kg) for any reason.
- Always assist the sick person during handling operations.
- Completely lock the chair wheels each time the patient is transferred.
- Before transferring the patient to/from the chair, make sure that the chair is placed on a flat, stable surface with the wheels locked.
- Before leaning forward or getting up from the chair, make sure that the front wheels are straight and the rear wheels are locked, to prevent a possible tip-over.
- To avoid possible tipping, do not lean over the back of the chair.
- Keep the chair away from direct heat sources.

USE

Getting on the shower chair

- Bring the wheelchair as close as possible so that it is positioned next to the patient.
- Lock the rear wheel brakes.
- Lean against the frame structure by forcing yourself on your arms, then slowly descend towards the seat.

Getting off the shower chair

- Get as close as possible to a new seat.
- Lock the brakes of the rear wheels;
- Place your arms on the frame of the seat on which you are to move and slowly slide towards the seat as you stand up.

⚠ ATTENTION!

Assistance or supervision by experienced and qualified personnel is required to prevent slipping and falls.

WHEELS

The 4 wheels of the shower chair are adjustable in height by moving the clips into the holes on them.

BRAKES

The brakes are applied by pressing the lever located on the wheels of the chair. When adjusted correctly, the brake lever must fit completely against the wheel.

TOILET BASIN

Use the toilet bowl with the chair fully braked.

ADJUSTING THE TILT OF THE SHOWER CHAIR

To tilt the chair back, pull the lever installed on the push handle and push down.

To return the chair to its normal position, simply pull the lever installed on the handle of thrust. The hydraulic pistons are designed to return to the starting position without the user applies no force.

FOOTPEG HEIGHT ADJUSTMENT

The footrests of the chair can be adjusted to the desired height:

- Loosen the adjustment knob;
- Adjust the footplate to the necessary height;
- Tighten the knob securely. Always check that the two footrests are at the same height.

MAINTENANCE

The shower chair requires a few elementary periodic maintenance, below we list the checks to be carried out periodically.

1. Check the brakes: monthly.
2. Cleaning of parts (seat, backrest, armrests, wheels, etc.) as needed.
3. Pin lubrications: monthly.

It is of fundamental importance that the shower chair is completely checked in all its parts when used for the first time, in particular that the screws and nuts are all tightened.

CLEANING AND DISINFECTION

Cleanliness

Use a cloth soaked in warm water and mild soap and dry with a clean, dry cloth.

Disinfection

If necessary, disinfect the product using a mild disinfectant detergent for normal household use. The toilet bowl should be regularly cleaned and disinfected using a normal bath disinfectant.

ATTENTION!

Do not use abrasives, acids, alcohol, chlorine-based detergents, disinfectants and acetone, as these substances cause abrasion of the shower chair components.

The manufacturer will not be responsible for damage caused by the use of materials that may deteriorate the surface of the product or by chemical corrosives used during cleaning.

GENERAL DISPOSAL CONDITIONS

When disposing of the device, never use normal municipal solid waste delivery systems. It is recommended to dispose of the device through the common municipal ecological islands for the planned recycling operations of the materials used.

WARRANTY

All Demarta products are warranted against defects in material or workmanship for a period of 2 (two) years from the date of sale of the product, subject to any exclusions and limitations specified below. This warranty is invalid in the event of misuse, abuse or modification of the product and for failure to follow the operating instructions. Demarta does not warrant the products against damage or defects in the following conditions: natural disasters, failure to follow the guidelines and instructions for use, or failure to perform maintenance as indicated in the manual.

Components that can wear and tear are not covered by this warranty if the damage is due to normal use of the product.