

# SUCTION/GAS/AIR FILTER

BACTERIAL/VIRAL PERFORMANCE

CASING: TRANSLUCENT STYROLUX

99.9999% BACTERIAL AND VIRAL  
EFFICIENCY NELSON LABORATORIES USA

EFFECTIVE FILTRATION AREA = 24.6 cm<sup>2</sup>

MAXIMUM OPERATING TEMPERATURE = 121 DEGREES C

MAXIMUM OPERATING PRESSURE = 20 psi

FILTERS DOWN TO .027 MICRON

FOR STERILE APPLICATIONS ETO  
(ETHYLENE OXIDE) IS STANDARD

DATE OF ISSUE 18/04/01

