

Nome Prodotto

Carrozzina da Transito Viky Transit Plus.

Destinazione d'uso

Le carrozzine sono un mezzo di trasporto indispensabile per le persone disabili con assente o ridotta capacità di deambulazione.

Lo scopo è quello di facilitare gli spostamenti e di salvaguardare il benessere fisico della persona, evitando la comparsa di complicanze correlate alla posizione seduta obbligatoria (piaghe da decubito, dolori posturali cronici...).



Codice	Taglia	Lunghezza	Larghezza	Altezza	Larghezza Piegata	Larghezza Seduta	Profondità Seduta	Altezza Schienale	Altezza Seduta da Terra	Portata Massima (Kg)	Dimensioni Ruote Anteriori	Dimensioni Ruote Posteriori	Peso Netto con Appoggiatesta	Dimensioni Scatola
Misure in Cm														
01.04.320D41FD	41	69	53.5	92	22.5	41	40	44	50	120	20 (8")	30 (12")	15	73.5x25.5x93
01.04.320D43FD	43	69	55.5	92	23	43	40	44	50	120	20 (8")	30 (12")	15.1	73.5x25.5x93
01.04.320D46FD	46	69	57.5	92	23.5	46	40	44	50	120	20 (8")	30 (12")	15.2	73.5x25.5x93
01.04.320D48FD	48	69	59.5	92	24	48	40	44	50	120	20 (8")	30 (12")	15.3	73.5x25.5x93

Caratteristiche

Telaio	Acciaio Verniciato a Polvere
Crociera	Doppia
Colore	Grigio Canna di Fucile (Opaco)
Ruote Anteriori	20 x 2,5 cm (8" x 1") PVC Grigie
Ruote Posteriori	30 cm (12") PU Grigie
Braccioli	Ribaltabili, Modello a Scrivania con Poggiatesta in PVC
Pedane	Estraibili e Regolabili con Supporto per Gambe
Poggiatesta	PVC Nero
Sedile e Schienale	Nylon Morbido Imbottito, lavabile e disinfettabile

Manutenzione e Pulizia

Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica

Conservazione e Smaltimento

Si consiglia l'uso ad un solo utente

Conservare al riparo da fiamme libere e fonti di calore

Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata

Confezione e Nr. Articoli

Scatola	1 Pezzo
---------	---------

Garanzia

Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
---------	--

Riferimenti del Prodotto

Fabbricante	Demarta Virgino S.a.S.
Classe di Dispositivo 93/42	Classe I
Codice ISO	12.22.18.003
Codice CND	Y122106
Codice RDM	1563892/R

Product Name















Transport Wheelchairs Viky Transit Plus.

Use

The wheelchairs are an essential way of transport for disabled people with none or limited walking ability.

The purpose is to facilitate movements and to safeguard the physical well-being of the person, avoiding the appearance of complications related to the sitting position (bed sores, chronic postural pains ...).



Code	Size	Length	Width	Height	Width Folded	Seat Width	Seat Depth	Back Height	Floor to Seat Height	Weight Capacity	Front Wheel	Rear Wheel	Weight with Footrest	Carton Size
Measures in Cm														
01.04.320D41FD	41	69	53.5	92	22.5	41	40	44	50	120	20 (8")	30 (12")	15	73.5x25.5x93
01.04.320D43FD	43	69	55.5	92	23	43	40	44	50	120	20 (8")	30 (12")	15.1	73.5x25.5x93
01.04.320D46FD	46	69	57.5	92	23.5	46	40	44	50	120	20 (8")	30 (12")	15.2	73.5x25.5x93
01.04.320D48FD	48	69	59.5	92	24	48	40	44	50	120	20 (8")	30 (12")	15.3	73.5x25.5x93

Features

Frame	Powder Coating Steel Frame
Cross Bar	Double Cross Brace
Color	Smokey Grey (Dull Polish)
Front Wheels	20 x 2,5 cm (8" x 1") PVC Grey Tyre
Rear Wheels	30 cm (12"), Spoke Wheels, PU Grey Tyre
Armrest	Flip Back Desk Model Armrest with PVC Armpad
Footrest	Swing Away and Length Adjustable with Leg Support
Footplate	Plastic Black
Seat and Backrest	Soft and Padded Nylon, washable and disinfected

Maintenance and Cleaning

Disinfection is possible with 5% diluted disinfectants even with alcohol base

Storage and Disposal

Single-user Product

Keep away from heat, sparks, open flame or other sources of ignition

Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

Package and Nr. Items

Box 1 Piece

Warranty

European 24 months for manufacturing defects with the exception of consumable parts

Product References

Manufacturer	Demarta Virgino S.a.S.
Directive 93/42	Class I
ISO Code	12.22.18.003
CND Code	Y122106
RDM Code	1563892/R