

**Nome Prodotto**

Protesi Maxima Simmetrica Piatta.



| Codice         | Taglia | Taglia Francese | Taglia Europea |
|----------------|--------|-----------------|----------------|
| 06.01.4001P.0  | 0      | 65              | 50             |
| 06.01.4001P.1  | 1      | 70              | 55             |
| 06.01.4001P.2  | 2      | 75              | 60             |
| 06.01.4001P.3  | 3      | 80              | 65             |
| 06.01.4001P.4  | 4      | 85              | 70             |
| 06.01.4001P.5  | 5      | 90              | 75             |
| 06.01.4001P.6  | 6      | 95              | 80             |
| 06.01.4001P.7  | 7      | 100             | 85             |
| 06.01.4001P.8  | 8      | 105             | 90             |
| 06.01.4001P.9  | 9      | 110             | 95             |
| 06.01.4001P.10 | 10     | 115             | 100            |
| 06.01.4001P.11 | 11     | 120             | 105            |
| 06.01.4001P.12 | 12     | 125             | 110            |
| 06.01.4001P.13 | 13     | 130             | 115            |
| 06.01.4001P.14 | 14     | 135             | 120            |

**Destinazione d'uso**

Le protesi mammarie esterne vengono utilizzate in soggetti che hanno subito una mastectomia. Grazie a questi dispositivi la morfologia del corpo resta invariata.

**Descrizione**

La PROTESI MAXIMA SIMMETRICA PIATTA è composta da un gel di silicone avvolto da un sottile film di poliuretano; concepita per ripristinare una visione anatomico-estetica del corpo. La protesi è anallergica ed esteticamente naturale, la forma a goccia la rende adatta ad ogni tipo di seno.

| <b>Materiali</b>   |  |
|--|--|
| Protesi  | Silicone   |
| <b>Manutenzione e Pulizia</b>  |  |
| La pulizia è possibile con un panno umido e un detergente neutro   |  |
| Si consiglia di utilizzare deodoranti, profumi o altri cosmetici prima di indossare le protesi                       |  |
| <b>Conservazione e Smaltimento</b>   |  |
| Si consiglia l'uso ad un solo utente   |  |
| Conservare al riparo da fiamme libere e fonti di calore  |  |
| Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata |  |
| <b>Confezione e Nr. Articoli</b>   |  |
| Scatola  | 1 Pezzo  |
| <b>Garanzia</b>  |  |
| Europea  | 24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili |
| <b>Riferimenti del Prodotto</b>  |  |
| Fabbricante  | Demarta Virginio S.a.S.  |
| Classe di Dispositivo 93/42  | Classe I   |
| Codice ISO   | 06.30.18.006   |
| Codice CND   | Y063018  |
| Codice RDM   | 1463825/R  |

**Product Name**

Symmetrical Prosthesis Maxima with Flat Back.



| Code           | Size | French Size | European Size |
|----------------|------|-------------|---------------|
| 06.01.4001P.0  | 0    | 65          | 50            |
| 06.01.4001P.1  | 1    | 70          | 55            |
| 06.01.4001P.2  | 2    | 75          | 60            |
| 06.01.4001P.3  | 3    | 80          | 65            |
| 06.01.4001P.4  | 4    | 85          | 70            |
| 06.01.4001P.5  | 5    | 90          | 75            |
| 06.01.4001P.6  | 6    | 95          | 80            |
| 06.01.4001P.7  | 7    | 100         | 85            |
| 06.01.4001P.8  | 8    | 105         | 90            |
| 06.01.4001P.9  | 9    | 110         | 95            |
| 06.01.4001P.10 | 10   | 115         | 100           |
| 06.01.4001P.11 | 11   | 120         | 105           |
| 06.01.4001P.12 | 12   | 125         | 110           |
| 06.01.4001P.13 | 13   | 130         | 115           |
| 06.01.4001P.14 | 14   | 135         | 120           |

**Use**

External breast prosthesis are used in subjects who have undergone mastectomy. Thanks to these devices the morphology of the body remains unchanged.

**Description**

The SYMMETRICAL PROSTHESIS MAXIMA WITH FLAT BACK is composed of a silicone gel wrapped in a thin polyurethane film; designed to restore an anatomical-aesthetic vision of the body. The prosthesis is hypoallergenic and aesthetically natural, the teardrop shape makes it suitable for all types of breasts.

| Materials  |  |
|--|--|
| Prosthesis   | Silicone Gel   |
| Maintenance and Cleaning   |  |
| Cleaning is possible with a damp cloth and a mild soap   |  |
| Avoid direct contact with cosmetics (deodorant, perfumes...), please use them before putting the prosthesis on         |  |
| Storage and Disposal   |  |
| Single-user Product  |  |
| Keep away from heat, sparks, open flame or other sources of ignition   |  |
| Dispose of the device in compliance with current regulations on environmental protection and separate waste collection |  |
| Package and Nr.Items   |  |
| Box  | 1 Piece  |
| Warranty   |  |
| European   | 24 months for manufacturing defects with the exception of consumable parts |
| Product References   |  |
| Manufacturer   | Demarta Virginio S.a.S.  |
| Directive 93/42  | Class I  |
| ISO Code   | 06.30.18.006   |
| CND Code   | Y063018  |
| RDM Code   | 1463825/R  |