

**Nome Prodotto**

Linea di Contenzione ADIGGY: Supporto per Abduzione.

Codice	Taglia
01.12.ADX30	Unica





**Destinazione d'uso**

La Linea di Contenzione ADIGGY comprende tutti quei ausili utilizzati per immobilizzare il paziente alla sedia a rotelle, letto ecc. Vengono impiegati per garantire la massima sicurezza di chi li utilizza per evitare cadute o traumi di varia natura.

**Descrizione**

Il SUPPORTO PER ABDUZIONE mantiene gli arti inferiori in una posizione leggermente divaricata e, inoltre, evita lo scivolamento dalla sedia a rotelle.

Il fissaggio avviene tramite fibbie con chiusura a clips.

<b>Materiali</b>	
Supporto per Abduzione	Nylon
<b>Manutenzione e Pulizia</b>	
	Lavaggio Max 30° C
	Non Stirare
	Non Candeggiare
	Non Utilizzare l'Asciugatrice
<b>Conservazione e Smaltimento</b>	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
<b>Confezione e Nr. Articoli</b>	
Scatola	1 Pezzo
<b>Garanzia</b>	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
<b>Riferimenti del Prodotto</b>	
Fabbricante	Demarta Virginio S.a.S.
Classe di Dispositivo 93/42	Classe I
Codice CND	V0899
Codice RDM	204439/R



**Product Name**

ADIGGY Line: Abduction Support.

Code	Size
01.12.ADX30	One





**Use**

The ADIGGY Line includes all the aids used to immobilize the patient in a wheelchair, bed, etc. They are used to ensure maximum safety of those who use them to prevent falls or traumas of various kinds.



**Description**

The ABDUCTION SUPPORT keeps the lower limbs in a slightly apart position and prevents slipping from the wheelchair. The fixing is done with plastic buckles.

Materials	
Abduction Support	Nylon
Maintenance and Cleaning	
	Water temperature not above 30° C
	Do not iron
	Do not bleach
	Do not dry in dryer
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.a.S.
Directive 93/42	Class I
CND Code	V0899
RDM Code	204439/R